

Ultrasound Reports

Staff and Doctor Guide

Open the dashboard	Double-click Ultrasound Dashboard on the Desktop (or open http://localhost:8000 in the browser)
Something stopped working	Double-click Restart Ultrasound Services on the Desktop
Print a report	Open the patient, click PDF / Print

1. How it works

After a scan, the ultrasound machine sends the images to this PC. About half a minute after the last image, the system makes the patient's report by itself: the correct report form, with the patient's details, the measurements from the scanner and the scan images already filled in.

The doctor then opens the report on the dashboard, checks it, completes the findings and impression, saves it and prints it. Nothing has to be copied or typed again.

Tip: the report is ready only once the scanner has **sent** the images. If a patient is missing, check that the scan was sent from the ultrasound machine.

2. The day's list

The dashboard opens on **today's scans**, the newest patient at the top. Each patient has a number for the day (1 = first patient). The list updates by itself when a new report is ready, and moves on to the new day after midnight.

The status of each report

Status	Meaning
Awaiting doctor	Report made, not yet checked and saved by the doctor.
Reviewed (green)	The doctor has edited and saved the report.
Orange warning line under the name	Needs checking: something needs attention, and the line says what (see section 7).
Measurements only	The type of scan could not be recognised, so no Word report was made. Choose the report type on the patient's page (section 5).

Finding a patient

- **Other days:** use **Previous** / **Next** at the top, the date box, or the calendar on the left (year, month, day).
- **Search:** type a name or patient ID in the search box at the top left.
- **Open a patient:** click anywhere on the row.

3. Working on a report

The patient's page shows the report as it will print, on the left, and the scan images with their measurements on the right.

Typing in the report

- 1 Click in the report and type, exactly as in Word. **B**, *I* and U on the toolbar make text bold, italic or underlined.
- 2 A value the scanner did not send is printed as _____. Fill in or delete every _____ before printing.
- 3 Click **Save** (or press **Ctrl+S**). The status changes to **Reviewed**.

Tip: the grey lines marked **Page 2**, **Page 3** show where a new printed page will start. Unsaved changes are shown as **Unsaved changes** next to the Save button.

Putting values from the scan into the report

- 1 Click in the report where the value should go (for example on the ____ after *Liver*).
 - 2 In the panel on the right, click the value under the matching image. It is typed into the report.
- **Suggested from the scan** (abdomen and pelvis scans): the system suggests which organ each measurement belongs to. Click a suggestion to see its image and check it, then click **Fill**, or **Fill all blanks**.
 - **Calculate CPR** (growth scans): after filling the umbilical artery and MCA PI, this button works out the CPR and puts it in the report.

Moving between patients

Use **Previous** / **Next** at the top of the page, or **Alt + Up arrow** / **Alt + Down arrow** on the keyboard.

4. Images in the report

Every image of the scan is shown on the right. The report's last page (*Ultrasound images*) holds the images chosen for print - the first six taken, to start with.

- **+ Add to report** under an image puts it on the report's images page.
- Images already in the report have a **green border** and the button **In report - remove** (with a tick), which takes them out.
- The images page is rebuilt straight away, in the order the images were taken.

5. Wrong report type

If the wrong exam was selected on the scanner (for example an abdomen scan made on the obstetric setting), the report is made on the wrong form. To correct it:

- 1 On the patient's page, choose the right form in **Report type** (top right).
- 2 Click **Remake report** and confirm. The report is made again on the new form.

Note: changes already typed into the old report are **not** carried over. Remake first, then edit. Remaking is possible until about 30 days after the scan; once the report has been reviewed, the original scan is then cleared from the system to save space.

6. Printing

- 1 Open the patient and click **PDF / Print**. The button shows **Preparing...** for a few seconds.
 - 2 The print window opens over the page. Choose the printer and click **Print**.
 - 3 The window closes by itself; you stay on the patient's page.
- **To keep a PDF** (for example to send by WhatsApp or e-mail): in the print window choose **Save as PDF** as the printer.
 - **Word file** (download arrow) gives the report as a Word document.
 - Reports print on the **pre-printed letterhead**: load letterhead paper in the printer.
 - Unsaved changes are saved automatically before printing.

Tip: if page 2 comes out before page 1, the printer is set to print the last page first. Change **Page order** in the printer's settings (Windows Settings > Printers > Printing preferences).

7. Warnings (Needs checking)

A report with an orange warning line is complete but has something to check. The line says what, on the list and on the patient's page.

Warning	What to do
scan dated ... but received ...	The ultrasound machine's date is wrong (often after a power cut). Correct the date and time on the machine. The report is filed under the wrong day.
no measurements extracted / measurement(s) left blank	Fill the values in by hand from the images on the right.
scan type not identified / low classification confidence	Check the report form is right; if not, use Remake report (section 5).
doctor name not set in Settings	Enter the doctor's name in Settings (section 8). Until then the declaration line is blank.
report folder ... not reachable	The folder for saving reports could not be used; the report was saved on this PC instead. Check the drive or folder in Settings .

8. Settings

Click **Settings** in the dashboard.

- **Reporting doctor** - name and qualification, printed in the PCPNDT declaration and the signature on every report.
- **Default "Ref by"** - used when no referring doctor was entered on the scanner.
- **Report folder** - where the Word and PDF reports are saved (for example D:\GeneratedReports\reports). Leave it blank to keep them on this PC.

Changing the report forms

Each report type has a Word form in **C:\scanning-machine\data\forms** (one file per type, for example *growth_scan.docx*). Open it in Word to change the standard wording; new reports use the changed form.

Careful: do not change or remove the words in double curly brackets, such as **{{bpd}}** or **{{patient_name}}** - that is where the system puts the values. Keep a copy of a form before changing it.

9. For the sonographer

The more the scanner sends, the less has to be typed:

- Select the **correct exam / preset** for each patient (obstetric, abdomen, gynae, fertility, breast).
- Enter the **patient's name and age** (for example *LAKSHMI 27*) and, for obstetric scans, the **LMP**.
- Use the **labelled measurement tools** (ovary, uterus, follicle, BPD...) rather than the plain distance calliper - labelled values go straight into the report.
- Keep the machine's **date and time** correct.
- **Send** the study to the PC at the end of each scan.

10. When something goes wrong

Problem	What to do
The dashboard does not open, or no new reports appear	Double-click Restart Ultrasound Services on the Desktop and click Yes . Wait about a minute. It shows OK or FAILED for each part and opens the dashboard when all are OK.
A patient's report is missing	Check the scan was sent from the ultrasound machine. Wait a minute, then check the list again.
PDF / Print says <i>Could not create the PDF</i>	Close Word if the report is open in it and try again. If it still fails, run C:\scanning-machine\deploy\dashboard_as_user.bat and press Enter, or click Word file and print from Word.
The report is on the wrong day	The ultrasound machine's date is wrong - correct it (see section 7).
Save says <i>The report is open in Word</i>	Close the report in Word, then save again.
Restart shows FAILED	Restart the PC. If it still fails, call the person who installed the system and tell them which part failed.

Where things are kept

- **Reports** (Word and PDF) and the scan images: in the report folder set in Settings.
- The original scan images stay on the PC for about **30 days**; once the report is reviewed they are then cleared to save space. The report and the images shown on the dashboard are kept.
- Keep the PC switched on and connected during clinic hours: a sleeping or switched-off PC receives no scans.

Patient reports are confidential medical records (PCPNDT). Do not share them outside the clinic except with the patient and the referring doctor.